



Revolutionizing
Behavior Research



AnnoStarMobile

Features:

- Manual Behavior Annotation System for Tablets and Mobile Phones!
- Currently available for **Android 2.2, 3.0, 3.1 & 3.2**
- Annotation based on touchscreen button presses
- Multiple types of button presses are possible to detect different types of events such as instantaneous, fixed-range, and variable-range events
- Real-time or offline scoring
- Scoring with or without video input
- Full Integratable with our Clever Sys products that allows both Automated AND Manual scoring to be performed simultaneously
- Detailed statistics about the manual annotated events that occurred during the experiment

Capabilities of AnnoStar:

- Touchscreen Button based Annotation
- Instantaneous Events
- Fixed-range events
- Variable-range events
- Manual events are handled just like automated events allowing full statistical analysis.

Requirements:

- Tablet running Android Version 2.2, 3.0, 3.1 or 3.2
- Mobile Phone running Android 2.2

AnnoStarMobile is our newest addition to our product family and is a mobile version of our versatile **AnnoStar** software. AnnoStarMobile has now available for the **Android platform** (tablets and mobile phones). It is an extremely user-friendly, rapid behavior annotation software. It allows the user to view the video of animals on the screen and manually annotate their behaviors along side, with ease. This Mobile App contains touchscreen buttons for rapidly annotating behaviors.



The software contains the following convenient features:

- Using touchscreen buttons as a mechanism for registering events.
- Ability to define a custom set of behaviors as annotation buttons. This allows the software to be used for any imaginable annotation task!
- Ability to export annotations to statistics software.
- Ability to operate with or without a video feed.
- Without a video feed, the App can be used to simply log tasks the user is scoring live. For example, the user can log events of zoo animals while walking around the zoo.
- With a saved video file, the user can load the video into the App, and playback the video on the mobile device. During this playback, the user can annotate the events. The allows the user to be performing annotations of stored video virtually anywhere – on the beach, while commuting on a train, at home, etc!!
- With an optional external camera connected to the device, one can even annotate using live-video feed!

Currently, the Android App is available for purchase either standalone for a mobile device or as an add-on option to our full PC-based AnnoStar software. Apps cannot be shared across mobile devices. If you have multiple mobile devices that you would like to use, you will need to purchase an App for each device. A complete turn-key system including a compatible Tablet/Phone and the App are also available for purchase from Clever Sys.